

DNA extraction from Tissue (mouse tail)

Preparation of test samples

Cut 3~5mm of the mouse tail (c.a.10mg) and put it in the micro test tube (usually use 0.2ml or 0.5ml tubes for PCR)



Reagents

Mix the CellEase A, B and distilled water (20 µl CellEase A, 20 µl CellEase B, 60 µl distilled water)



Add 100 µl of the mixture to the samples.



Incubate at 72°C for 6 minutes
Then incubate at 94°C for 3 minutes



Transfer 6 µl of extracts to PCR reaction mixture and amplify the target DNA fragment

PCR

6.0 µl	Test sample
5.0 µl	x 10 buffer (+Mg ²⁺)
5.0 µl	dNTPs
1.0 µl	Forward Primer (10pmol/µl)
1.0 µl	Reverse Primer (10pmol/µl)
0.5 µl	Ex Taq (5 U/µl)

Fill up to 50 µl by distilled water


PCR Cycle

94°C	1min
94°C	30sec
55°C	30sec
72°C	60sec
72°C	4min

35 Cycles

Primers : Mouse β-globin gene (494bp)

<Comparison with conventional CellEase kit >



M	Marker (100bp ladder)	
1	DNA extract by using CellEase Tissue II	conc.
2	DNA extract by using CellEase Tissue II	× 10
3	DNA extract by using CellEase Tissue II	× 100
4	DNA extract by using conventional CellEase	conc.
5	DNA extract by using conventional CellEase	× 10
6	DNA extract by using conventional CellEase	× 100

- ※ The protocol of conventional CellEase and CellEase II kit were followed by the instruction manual respectively.
- ※ The thickness of the DNA bands were depending upon the amount of test samples and a parts of tissue including bone, skin or fat.

The clear DNA bands were detected from more than × 10 dilution of DNA extracts by using CellEase Tissue II